

NAD-Malate Dehydrogenase (NAD-MDH) Activity Assay Kit

Note: It is necessary to predict 2-3 large difference samples before the formal determination.

Operation Equipment: Spectrophotometer

Cat No: AK0498

Size: 50T/48S

Components:

Extract solution: 50 mL×1, store at 4C .

Reagent I: 40 mL×1, store at 4C .

Reagent II: Powder×2, store at -20C . Add 360 μL of distilled water when the solution will be used. The rest of reagentstore at -20C .

Reagent III: Powder×2, store at -20C . Add 327 μL of distilled water when the solution will be used. The rest of reagentstore at -20C .

Product Description:

Malate Dehydrogenase (MDH, EC 1. 1 . 1 .37) is widely exist in animal, plant, microbe and cells culture. MDH in mitochondria is one of the key enzymes of TCA cycle, which catalyzes the formation of oxaloacetic acid from malic acid. In contrast, MDH in serum catalyzes the formation of malic acid from oxaloacetic acid. Oxaloacetic acid is an important intermediate that connects several important metabolic pathways. MDH plays an important role in many physiological activities of cells, including mitochondrial energy metabolism, malic acid-aspartic acid shuttle system, reactive oxygen species metabolism and disease resistance. According to the different coenzyme specificity, MDH is divided into NAD- dependent MDH and NADP- dependent MDH. Bacteria usually only contain NAD-MDH. NAD-MDH is distributed in the cytoplasm and mitochondria in eukaryotic cells.

NAD- MDH catalyzes NADH to reduce oxaloacetic acid into malic acid, which resulting in a decrease in absorption at 340nm.

Required but Not Provided:

Ultraviolet spectrophotometer, desk centrifuge, water-bath, adjustablepipette, 1 mL quartzcuvette, distilled water.

Protocol

I. Preparation:

1. Cells or bacterial

Collect bacteria or cells into the centrifuge tube. The liquid in the upper layer is discarded after centrifugation. It is suggested to take about 2 million bacteria/cell and add 400 μ L of Extract solution. Bacteria/cell is split by ultrasonic (power 20%, ultrasonic 3s, interval 10s, repeat for 30 times). Centrifuge at 8000 \times g and 4°C for 10 minutes. Take the supernatant and place it on ice for testing.

2. Tissue:

It is suggested to take about 0.05 g of tissue and add 1 mL of Extract solution. Then ice bath homogenization should be carried out. Centrifuge at 8000 × g and 4°C for 10 minutes. The supernatant is placed on ice for testing.

3. Serum: detect directly.

II. Procedure:

Preheat ultraviolet spectrophotometer for 30 minutes, adjust wavelength to 340 nm, set zero with Distilled water. Preheat Reagent I in 37°C for 15 minutes.

Reagent (μL)	Test tube (T)	Blank tube (B)
Sample	20	-
Distilled water	-	20
Reagent I	760	760
Reagent II	10	10
Reagent III	10	10

The above reagents are added into a 1 mL of quartz cuvette in sequence. After fully mixing, the initial absorbance A₁ and the absorbance A₂ after 1 minute of reaction are recorded at 340 nm wavelength, and the reaction temperature is kept at 37°C as far as possible. ΔA = A₁ - A₂. Record ΔA_T, ΔA_B.

III. NAD-MDH Calculation:

1. Serum

Unit definition: One unit of enzyme activity is defined as the amount of enzyme consumes 1 nmol of NADH per minute in the reaction system every milliliter serum.

$$\text{NAD-MDH (U/mL)} = (\Delta A_T - \Delta A_B) \div (\epsilon \times d) \times V_{RT} \div V_{SA} \div T = 6430 \times (\Delta A_T - \Delta A_B)$$

2. Tissue

1) Protein concentration:

Unit definition: One unit of enzyme activity is defined as the amount of enzyme consumes 1 nmol of NADH per minute in the reaction system every milligram tissue protein.

$$\text{NAD-MDH (U/mg prot)} = (\Delta A_T - \Delta A_B) \div (\epsilon \times d) \times V_{RT} \div (V_{SA} \times C_{pr}) \div T = 6430 \times (\Delta A_T - \Delta A_B) \div C_{pr}$$

2) Sample weight:

Unit definition: One unit of enzyme activity is defined as the amount of enzyme consumes 1 nmol of NADH per minute in the reaction system every gram tissue.

$$\text{NAD-MDH(U/g)} = (\Delta A_T - \Delta A_B) \div (\epsilon \times d) \times V_{RT} \div (W \div V_E \times V_{SA}) \div T = 6430 \times (\Delta A_T - \Delta A_B) \div W$$

3. Cells or germ

1) Protein concentration:

Unit definition: One unit of enzyme activity is defined as the amount of enzyme consumes 1 nmol of NADH per minute in the reaction system every milligram tissue protein.

$$\text{NAD-MDH(U/mg prot)} = (\Delta A_T - \Delta A_B) \div (\epsilon \times d) \times V_{RT} \div (V_{SA} \times C_{pr}) \div T = 6430 \times (\Delta A_T - \Delta A_B) \div C_{pr}$$

2) Cells or germ

Unit definition: One unit of enzyme activity is defined as the amount of enzyme consumes 1 nmol of

NADH per minute in the reaction system every 10 thousand cells or germ.

$$\text{NAD-MDH(U/ } 10^4 \text{ cell)} = (\Delta A_T - \Delta A_B) \div (\epsilon \times d) \times V_{RT} \div (500 \times V_{SA}) \div T = 13 \times (\Delta A_T - \Delta A_B)$$

V_{RT} : Total reaction volume, 0.8mL;

V_{SA} : Sample volume, 0.02mL;

ϵ : NADH molar extinction coefficient, 6.22×10^{-3} mL/nmol/cm;

d: Light path of cuvette, 1cm;

T: Reaction time, 1 minute;

Cpr: Protein concentration, mg/mL;

500: Cells or germ, 10^4 /mL

Note:

1. The extraction of crude enzyme must be done at 0C-4C to prevent enzyme denaturation and inactivation.
2. Reagent II, III and sample must keep on ice during experiment in order to avoid denaturation and inactivation.
3. When micro quartz cuvette with the light diameter of 1 cm is used, if the absorbance less than 0.7 or $\Delta A > 0.5$, dilute the liquid before recommending. When 96 well flat-bottom plate (UV plate) is used, if the absorbance less than 0.4 or $\Delta A > 0.3$, dilute the liquid before recommending.
4. It is suggested that two people to cooperate in the experiment, one adding the sample and the other colorimetric.
5. Blank tube only needs to measure 1-2 times.

Experimental instances

1. Take 0.05g of rat liver, add 1 mL of extract solution, homogenate and grind. Take the supernatant and detect according to the measured steps. Calculate $\Delta A_T = A_{1T} - A_{2T} = 0.917 - 0.293 = 0.624$, $\Delta A_B = A_{1B} - A_{2B} = 0.789 - 0.784 = 0.005$, calculate the enzyme activity according to sample weight:

$$\text{NAD-MDH (U/g weight)} = 6431 \times (\Delta A_T - \Delta A_B) \div W = 6431 \times 0.619 \div 0.05 = 79616 \text{ U/g weight.}$$

2. Take 0.05g Willow leaf, add 1 mL of extract solution, homogenate and grind. Take the supernatant and detect according to the measured steps. Calculate $\Delta A_T = A_{1T} - A_{2T} = 0.812 - 0.741 = 0.071$, $\Delta A_B = A_{1B} - A_{2B} = 0.789 - 0.784 = 0.005$, calculate the enzyme activity according to sample weight:

$$\text{NAD-MDH (U/g weight)} = 6431 \times (\Delta A_T - \Delta A_B) \div W = 6431 \times 0.066 \div 0.05 = 8489 \text{ U/g weight.}$$

3. Take 20 μ L serum of mouse to detect directly, calculate $\Delta A_T = A_{1T} - A_{2T} = 0.855 - 0.789 = 0.066$, $\Delta A_B = A_{1B} - A_{2B} = 0.789 - 0.784 = 0.005$, calculate the enzyme activity according to volume of serum:

$$\text{NAD-MDH (U/mL)} = 6431 \times (\Delta A_T - \Delta A_B) = 6431 \times 0.061 = 392 \text{ U/mL.}$$



References

[1] Yao Y X, Li M, Zhai H, et al. Isolation and characterization of an apple cytosolic malate dehydrogenase gene reveal its function in malate synthesis[J]. Journal of plant physiology, 2011, 168(5) : 474-480 .

Related products

AK0560/AK0559 Coenzyme I NAD (H) Content Assay Kit

AK0500/AK0499 NAD Kinase (NADK) Assay Kit

AK0528/AK0527 NADH oxidase (NOX) Activity Assay Kit

AK0484/AK0483 NAD Malic Enzyme (NAD-ME) Assay Kit