

Superoxide Dismutase (SOD) Activity Assay Kit

Note: It is necessary to predict 2-3 large difference samples before the formal determination.

Operation Equipment: Spectrophotometer

Cat No: AK0584

Size: 50T/24S

Components:

Extraction reagent: 60 mL×1. Storage at 4°C .

Reagent I : 15 mL×1. Storage at 4°C .

Reagent II : 160 μL×1. Storage at 4°C . Mix by pipetting after centrifugation.

Reagent III: 11 mL×1. Storage at 4°C .

Reagent IV: Powder×1. Storage at 4°C .

Reagent V: 2 mL×1. Storage at 4°C . Add Reagent IV to Reagent V before use and shake with an oscillator to mix thoroughly. It can be stored 3 months.

Product Description:

Superoxide dismutase (SOD, EC 1. 15. 1. 1) is widely found in animals, plants, microorganisms and cultured cells. It catalyzes the superoxide anion to form H₂O₂ and O₂ . SOD is not only the superoxide anion scavenging enzyme, but also the main H₂O₂ producing enzyme, which plays an important role in the biological antioxidant system.

Superoxide anion (O₂⁻) is produced by the xanthine and xanthine oxidase reaction system. O₂⁻ can reduce blue tetrazole to form blue formazan, which has absorbance in 560 nm. SOD can remove O₂⁻ and inhibit the formation of methionine. The darker the blue color of the reaction solution, the lower the activity of SOD. The lighter the blue color of the reaction solution, the higher the activity of SOD.

Reagents and Equipment Required but Not Provided:

Spectrophotometer, table centrifuge, transferpettor, 1 mL glass cuvette, mortar/homogenizer, ice and distilled water.

Operation steps:

I. Sample preparation:

1. Bacteria or cells: collecting bacteria or cells into the centrifuge tube, discard supernatant after centrifugation. According to the proportion of bacteria or cells (10⁴ cells): extraction solution volume (mL) of 1:5- 10 to extract. It is suggested that 5 million of bacteria or cells amount with 1 mL of Extraction reagent. Splitting the bacteria or cells with ultrasonication (placed on ice, ultrasonic power 200W or 20%, working time 3s, interval 10s, repeat for 30 times). Centrifuge at 8000 ×g for 10 minutes at 4°C to remove insoluble materials, and take the supernatant on ice before testing.
2. Tissue: according to the proportion of tissue weight (g): extraction solution volume (mL) of 1:5- 10 to extract. It is suggested that 0.1 g of tissue with 1 mL of Extraction reagent and fully homogenized on ice

3. bath. Centrifuge at 8000 ×g for 10 minutes at 4°C to remove insoluble materials, and take the supernatant on ice before testing.
4. Serum (plasma) sample: detect sample directly.

II. Determination procedure:

1. Preheat the spectrophotometer for 30 minutes, adjust wavelength to 560 nm and set zero with distilled water.
2. Keep Reagent I , Reagent III , Reagent V in water bath for more than 5 minutes at 37°C(mammal) or 25°C (other species).
3. Add reagents with the following list:

Reagent (μL)	Test tube (T)	Control tube (C)	Blank tube (B1)	Blank tube (B2)
Sample	90	90	-	-
Reagent I	240	240	240	240
Reagent II	6	-	6	-
Reagent III	180	180	180	180
Distilled water	480	486	570	576
Reagent V	30	30	30	30

Mix thoroughly and the mixture is incubated at room temperature for 30 minutes. Add the mixture into 1mL glass cuvette, and detect the absorbance value of each tube at 560 nm. $\Delta A_T = A_T - A_C$, $\Delta A_B = A_{B1} - A_{B2}$. If there is precipitation at the bottom, mix thoroughly and then measure.

III. Calculation:

1. Inhibition percentage:

$$\text{Inhibition percentage} = [\Delta A_B - \Delta A_T] \div \Delta A_B \times 100\%$$

The inhibition percentage should be in 30%~70% (the value close to 50% will have a more accurate result). If the calculated inhibition percentage is less than 30% or more than 70%, it is usually necessary to adjust the sample addition amount and re determine. If the percentage of inhibition is too high, the sample should be diluted properly. If the percentage of inhibition is too low, the sample should be reprepared with a higher concentration.

2. Unit definition: One unit of enzyme activity is defined as the amount of enzyme catalyzes the inhibition of 50% in the reaction system of the above xanthine oxidase.

3. Calculation

A. Serum (plasma) sample

$$\text{SOD (U/mL)} = [P \div (1-P) \times V_{rv}] \div V_s \times F = 11.4 \times P \div (1-P) \times F$$

B. Tissue, bacteria or cultured cells

- a) Protein concentration:

$$\text{SOD (U/mL prot)} = [P \div (1-P) \times V_{rv}] \div (V_s \times C_{pr}) \times F = 11.4 \times P \div (1-P) \div C_{pr} \times F$$

b) Sample weight

$$\text{SOD (U/g weight)} = [P \div (1-P) \times V_{rv}] \div (W \times V_s \div V_{sv}) \times F = 11.4 \times P \div (1-P) \div W \times F$$

c) Bacteria or cell amount

$$\text{SOD (U/10}^4 \text{ cell)} = [P \div (1-P) \times V_{rv}] \div (500 \times V_s \div V_{sv}) \times F = 0.0228 \times P \div (1-P) \times F$$

V_{rv} : Total reaction volume, 1.026 mL;

V_s : Sample volume, 0.09 mL;

V_{sv} : Extraction volume, 1 mL;

C_{pr} : Sample protein concentration, mg/mL;

W : Sample weight, g;

500: Total number of bacteria and cells, 5 million.

P : Inhibition percentage, %;

F : Sample dilution multiple.

Note:

1. The Sample and Reagent II should be placed on ice when using.
2. When there are many samples, the working solution (including Reagent I, II and III) can be configured according to the table. Reagent V must be added finally.
3. After the reaction completed, there may be precipitation formed, which can be determined after mixing.

Experimental Examples:

1. 0.1 g of Echinochloa crusgalli is added into 1 mL of Extraction reagent for homogenization. After the supernatant is taken, the operation is carried out according to the determination steps. The results showed that $\Delta A_T = A_T - A_C = 0.335 - 0.012 = 0.323$, $\Delta A_B = A_{B1} - A_{B2} = 0.957 - 0.003 = 0.954$. Inhibition percentage = $(\Delta A_B - \Delta A_T) \div \Delta A_B \times 100\% = 72\%$, and the enzyme activity is calculated according to the sample mass.

SOD activity (U/g mass) = $11.4 \times \text{Inhibition percentage} (1 - \text{Inhibition percentage}) \times W = 293.14$ U/g mass.

2. 1 mL of Extraction reagent is added to 0.1 g of rat spleen for homogenization. After the supernatant is taken, the operation is carried out according to the determination steps. The results showed that $\Delta A_T = A_T - A_C = 0.563 - 0.213 = 0.35$, $\Delta A_B = A_{B1} - A_{B2} = 0.957 - 0.003 = 0.954$, inhibition percentage = $(\Delta A_B - \Delta A_T) \div \Delta A_B \times 100\% = 63.31\%$

SOD activity (U/g mass) = $11.4 \times \text{Inhibition percentage} (1 - \text{Inhibition percentage}) \times W = 196.71$ U/g mass.

3. 10 million cells is extracted and centrifuged by adding 1 mL of Extraction reagent, and then the operation is performed according to the determination steps. The results are as follows: $\Delta A_T = A_T -$

$$A_C = 0.614 - 0.015 = 0.599, \Delta A_B = A_{B1} - A_{B2} = 0.944 - 0.005 = 0.939, \text{inhibition percentage} = (\Delta A_B - \Delta A_T) \times \Delta A_B \times 100\% = 36.21\%$$

$$\text{SOD activity (U/10}^4 \text{ cell)} = \text{Inhibition percentage} \div (1 - \text{Inhibition percentage}) \times V_{TS} \div (1000 \times V_S \div V_{TS}) = 0.0065 \text{ U/10}^4 \text{ cell.}$$

Recent Protect Citations:

[1] Beibei Li, Yang Ding, Xiuli Tang, et al. Effect of L-Arginine on Maintaining Storage Quality of the White Button Mushroom (*Agaricus bisporus*). *Food and Bioprocess Technology*. April 2019; 12: 563-574. (IF3.032)

[2] Wang Li, Shi Qinghai, Li Kai, et al. Oral administration of Ginkgolide B alleviates hypoxia-induced neuronal damage in rat hippocampus by inhibiting oxidative stress and apoptosis. *Iranian Journal of Basic Medical Sciences*. February 2019; 22(2): 140- 145. (IF1.854)

[3] Yanan Wang, Chengzhen Liang, Zhigang Meng, et al. Leveraging *Atriplex hortensis* choline monooxygenase to improve chilling tolerance in cotton. *Environmental and Experimental Botany*. June 2019; 162: 364-373. (IF3.712)

[4] Fangzhou Chen, Ying Zhao, Huizhao Chen. MicroRNA-98 reduces amyloid β -protein production and improves oxidative stress and mitochondrial dysfunction through the Notch signaling pathway via HEY3 in Alzheimer's disease mice. *International Journal of Molecular Medicine*. October 2018; 91- 102. (IF2.784)

[5] Yang Yang, Li Jing, Wei Cong, et al. Amelioration of nonalcoholic fatty liver disease by swertiamarin in fructose-fed mice. *Phytomedicine*. June 2019; 59. (IF4. 18)

References:

[1] Spitz D R, Oberley L W. An assay for superoxide dismutase activity in mammalian tissue homogenates[J]. *Analytical Biochemistry*, 1989, 179(1):8- 18.

[2] Masayasu M, Hiroshi Y. A simplified assay method of superoxide dismutase activity for clinical use[J]. *Clinica Chimica Acta*, 1979, 92(3):337-342.

Related Products:

AK0293/AK0292 Polyphenol Oxidase (PPO) Activity Assay Kit

AK0578/AK0577 Phenylalnine Ammonialyase (PAL) Activity Assay Kit

AK0580/AK0579 Catalase (CAT) Activity Assay Kit

AK0598/AK0597 Peroxidase (POD) Activity Assay Kit