

OVA Conjugated Hyaluronic Acid (HA) Instruction Manual

SDPA212Ge21

Pan-species (General)

FORMULA

USAGE

Reconstitute in PBS or others.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.

| | |
|--------------------------------|---|
| Source | Protein Conjugation |
| Original Molecular Mass | n/a |
| Buffer Formulation | PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% |

| | |
|--------------------------|---------------------------|
| Traits | Trehalose and Proclin300. |
| Purity | Freeze-dried powder |
| Isoelectric Point | > 95% |
| Applications | n/a |
| | Immunogen; SDS-PAGE; WB. |