

Eukaryotic Complement Component 5a (C5a) Instruction Manual

SFPA122Mu61

Mus musculus (Mouse)

Source	Eukaryotic expression
Host	293F Cell
Endotoxin Level	<1.0EU per 1 μ g (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	10.5kDa
Accurate Molecular Mass	14kDa(Analysis of differences refer to the manual)
Residues & Tags	Asn679~Arg755 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	9.3
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

NL HLLRQKIEEQ AAKYKHSVPK
KCCYD GARVN FYETCEERVA RVTIGPLCIR AFNECCTIAN KIRKESPHKP
VQLGR

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

