

Recombinant Interleukin 4 (IL4)

Instruction Manual

SIPA049Rb01

Oryctolagus cuniculus (Rabbit)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	15.2kDa
Accurate Molecular Mass	15kDa(Analysis of differences refer to the manual)
Residues & Tags	Arg26~Ser147 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.05% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 90%
Isoelectric Point	9.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
RGDII LPEVIKTLNI LTERKTPCTK  
LMIADALAVP KNTTTEREAVC RAATALRQFY LHHKVSWCFK EHGEGLDLRL  
LRGLDRNLCS MAKLSNCPGK EARQTTLEDF LDRLKTAMQE KYSKRQS
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

