

Recombinant Interleukin 6 (IL6)

Instruction Manual

SIPA051Si01

Rhesus monkey (Simian)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	24.0kDa
Accurate Molecular Mass	27kDa(Analysis of differences refer to the manual)
Residues & Tags	Pro27~Met212 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	7.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
PAPV LPGEDSKNVA APHSQPLTSS
ERIDKHIRYI LDGISALRKE TCNRSNMCES SKEALAENNL NLPKMAEKDG
CFQSGFNEDT CLVKIITGLL EFEVYLEYLQ NRFESSEEQA RAVQMSTKVL
IQFLQKKAKN LDAITTPPEPT TNASLLTKLQ AQNQWLQDMT THLILRSFKE
FLQSNLRALR QM
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

