

Recombinant Anti-Mullerian Hormone (AMH) Instruction Manual

SIPA104Hu03

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	61.1kDa
Accurate Molecular Mass	61kDa(Analysis of differences refer to the manual)
Residues & Tags	Leu19~Arg560 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	7.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
LL GTEALRAEEP AVGTSGLI FR EDLDWPPGSP
QEPLCLVALG GDSNGSSSPL RVVGALSAYE QAF LGAVQRA RWGPRDLATF
GVCNTGDRQA ALPSLRRLGA WLRDPGGQRL VVLHLEEV TW EPTPSLRFQE
PPPGGAGPPE LALLVLYPGP GPEVTVTRAG LPGAQSLCPS RDTRYLV LAV
DRPAGAWRGS GLALTLQPRG EDSRLSTARL QALLFGDDHR CFTRMTPALL
LLPRSEPAPL PAHGQLDTPV FPPRPSAEL EESPPSADPF LETLTRLVRA
LRVPPARASA PRLALDPDAL AGFPQGLVNL SDPAALERLL DGEEPLLLLLL
RPTAATTGDP APLHDPTSAP WATALARRVA AELQAAAAEL RSLPGLPPAT
APLLARLLAL CPGGPGGLGD PLRALLLLKA LQGLRVEWRG RDPRGPGRAQ
RSAGATAADG PCALRELSVD LRAERSVLIP ETYQANNCQG VCGWPQSDRN
PRYGNHVLL LKMQVRGAAL ARPPCCVPTA YAGKLLISLS EERISAHHVP
NMVATECGCR
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.

