

Recombinant Matrix Metalloproteinase 9 (MMP9) Instruction Manual

SIPA145Mu02

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	73.6kDa
Accurate Molecular Mass	74kDa(Analysis of differences refer to the manual)
Residues & Tags	Phe108~Ser508 with N-terminal His and GST Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 90%
Isoelectric Point	4.7
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

FQT FKGLKWDHHN ITYWIQNYSE DLPRDMIDDA FARAFVWGE
VAPLTFTRVY GPEADIVIQF GVAEHGDGYP FDGKDGLLAH AFPPGAGVQG
DAHFDDELW SLGKGVVIPT YYGNSNGAPC HFPFTFEGRS YSACTTDGRN
DGTPWCSTTA DYDKDGKFGF CPSELYTEH GNGEGKPCVF PFIFEGRSYS
ACTTKGRSDG YRWCATTANY DQDKLYGFCP TRVDATVVGG NSAGELCVFP
FVFLGKQYSS CTSDGRRDGR LWCATTSNFD TDKKWGFPCD QGYSFLVAA
HEFGHALGLD HSSVPEALMY PLYSYLEGFP LNKDDIDGIQ YLYGRGSKPD
PRPPATTTTE PQPTAPPTMC PTIPPTAYPT VGPTVGPTGA PSPGPTSSPS
PGPTGAPS

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.