

Recombinant Lipase, Hepatic (LIPC) Instruction Manual

SIPA268Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	16.4kDa
Accurate Molecular Mass	17/13kDa(Analysis of differences refer to the manual)
Residues & Tags	Ser347~Pro478 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	9.1
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
SPFK VYHYQFKIQF INQTETPIQT TFTMSLLGTK EKMQKIPITL  
KKGIASNKTY SFLITLDVDI GELIMIKFKW ENSAVWANVW DTVQTIIPWS  
TGPRHSGLVL KTIRVKAGET QQRMTFCSEN TDDLLLRP
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

TTTCCCTTCAGTTTATATATTCAGTAAAGCTGCTGTTATACCAATCTGACCCGATACAGTACTTTATATGTCATCTGGGAAAGAGCAATTCACAAATCCCTCATCTGEEAAGGCTTCTAGTAAAGAACTTCTCTATACGCTGGATGTGGTATCGGCGAGCTGATCATGTC
S P F K V Y B Y Q L K I Q F I N Q T E T P I Q T T F T W S I L G T K E R K Q K I P I T L G K G I A S N K T T S F L I T L D V D I G E L I N I :




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.