

Recombinant Cluster Of Differentiation 86 (CD86) Instruction Manual

SIPA354Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane
Predicted Molecular Mass	30.9kDa
Accurate Molecular Mass	26&31kDa(Analysis of differences refer to the manual)
Residues & Tags	Gln46~Cys282 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	8.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```

                                                                    QNQSL
SELVVFWDQD ENLVLNEVYL GKEKFDSVHS KYMGRTSFDS DSWTLRLHNL
QIKDKGLYQC IIHHKKPTGM IRIHQMNSL SVLANFSQPE IVPISNITEN
VYINLTCSSI HGYPEPKKMS VLLRTKNSTI EYDGMQKSQ DNVTELYDVS
ISLSVSFPDV TSNMTIFCIL ETDKTRLLSS PFSIELEDPQ PPPDHIPWIT
AVLPTVIICV MVFCLILWKW KKKKRPRNSY KC
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

