

Recombinant Glycogen Phosphorylase, Muscle (PYGM) Instruction Manual

SIPA356Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Cytoplasm
Predicted Molecular Mass	21.7kDa
Accurate Molecular Mass	26kDa(Analysis of differences refer to the manual)
Residues & Tags	Arg11~Gly187 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	5.6
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
                RKQISVRGLA  GVENVTELKK  NFNRLHFTL  VKDRNVATPR  
DYYFALAHTV  RDHLVGRWIR  TQQHYYEKDP  KRIYYLSLEF  YMGRTLQNTM  
VNLALENACD  EATYQLGLDM  EELEEIEEDA  GLGNGGLGRL  AACFLDSMAT  
LGLAAYGYGI  RYEFGIFNQK  ISGGWQMEEA  DDWLRYG
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

AGAAAGTAAATCAAGTGGTGGCTTACGGGATGGAAATGCTGGGACTTAAAGAACTTCAAGGACACTGCTTCACTGGTCAAGATGGCAATGCTATCGAGAGTACTATTTCGACTGGCCACTGTCCGGACACTGGTGGACCTGGATCCGACACAGAGACTTATCGAAAG
R K Q I S V R G L A G V E N V S E L K K N F N R E L H F T L V K D R N V A T P R D Y Y F A L A H T V R D E L V G R V I R T Q Q H Y Y E K I




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.