

Recombinant Ubiquitin (Ub) Instruction Manual

SIPA406Mi02

**Homo sapiens (Human), Mus musculus (Mouse), Rattus norvegicus (Rat),
Oryctolagus cuniculus (Rabbit), Canis familiaris; Canine (Dog), Sus scrofa; Porcine
(Pig), Bos taurus; Bovine (Cattle), Equus caballus; Equine (Horse)**

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus, Cytoplasm
Predicted Molecular Mass	29.5kDa
Accurate Molecular Mass	29kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Cys229 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	7.6
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MQIFVKLTG KTITLEVEPS DTIENVKAKI QDKEGIPPDQ QRLIFAGKQL  
EDGRTLSDYN IQKESTLHLV LRLRGGMQIF VKTLTGKTIT LEVPSDTIE  
NVKAKIQDKE GIPPDQQLI FAGKQLEDGR TLDYNIQKE STLHLVLRRLR  
GGMQIFVKTL TGKTITLEVE PSDTIENVKA KIQDKEGIPP DQQLIFAGK  
QLEDGRTLSD YNIQKESTLH LVLRLRGGC
```

USAGE

Reconstitute in ddH₂O to a concentration of 0.1-0.3 mg/mL. Do not vortex.


STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.