

Recombinant Thrombomodulin (TM) Instruction Manual

SIPA588Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	18.6kDa
Accurate Molecular Mass	25kDa(Analysis of differences refer to the manual)
Residues & Tags	Asp325~Asp481 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	4.3
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
DVDDCI LEPSPCPQRC VNTQGGFECH  
CYPNYDLVDG ECVEPVDPCF RANCEYQCQP LNQTSYLCVC AEGFAPIPHE  
PHRCQMFCNQ TACPADCNP TQASCECEP YILDDGFICT DIDECENGGF  
CSGVCHNLPG TFECICGPDS ALARHIGTDC D
```

USAGE

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

