

Recombinant 5T4 Instruction Manual

SIPA617Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane
Predicted Molecular Mass	38.8kDa
Accurate Molecular Mass	39kDa(Analysis of differences refer to the manual)
Residues & Tags	Pro34~Gln356 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	6.1
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
                P S S S V P S   S S T S P A A F L A
SGSAQPPPAE  RCPAACECE  AARTVKCVNR  NLLEVPADLP  PYVRNLFLTG
NQMTVLPAGA  FARQPPLADL  EALNLSGNHL  KEVCAGAFEH  LPGLRRLDLS
HNPLTNLSAF  AFAGSNASVS  APSPLEELIL  NHIVPPEDQR  QNGSFEGMVA
FEGMVAAALR  SGLALRGLTR  LELASNHFLF  LPRDLLAQLP  SLRYLDLRNN
SLVSLTYASF  RNLTHLESLH  LEDNALKVLH  NSTLAEWHGL  AHVKVFLDNN
PWVDCY MAD  MVAWLKETE V  VPDKARLTCA  FPEKMRNRGL  LDLNSSDLDC
DAVLPQ
```

USAGE

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.


STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.