

Recombinant Cluster Of Differentiation 28 (CD28) Instruction Manual

SIPA621Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane
Predicted Molecular Mass	18.7kDa
Accurate Molecular Mass	19kDa(Analysis of differences refer to the manual)
Residues & Tags	Asn20~Leu150 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	6.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
      N KILVKQSPLL VVDSNEVSL S CRYSYNLLAK  
EFRASLYKGV NSDVEVCVGN GNFTYQPQFR SNAEFNCDGD FNETVTFR L  
WNLHVNHTDI YFCKIEFMYP PPYLDNERSN GTIIHIKEKH LCHTQSSPKL
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

