

# Recombinant Ribonuclease H (RNASEH) Instruction Manual

SIPA704Hu01

**Homo sapiens (Human)**

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Cytoplasm
<b>Predicted Molecular Mass</b>	62.1kDa
<b>Accurate Molecular Mass</b>	62kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Met1~Asp286 with N-terminal His and GST Tag
<b>Buffer Formulation</b>	20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 97%
<b>Isoelectric Point</b>	9.3
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

## SEQUENCE

```
MSWLLFLAHR VALAALPCRR GSRGFGMFYA VRRGRKTGVF LTWNECRAQV  
DRFPAARFKK FATEDEAWAF VRKSASPEVS EGHENQHGQE SEAKASKRLR  
EPLDGDGHES AEPYAKHMKP SVEPAPPVSR DTFSYMGDFV VVYTDGCCSS  
NGRRRPRAGI GYVWGPQHPL NVGIRLPGRQ TNQRAEIHAA CKAIEQAKTQ  
NINKLVLYTD SMFTINGITN WVQGWKKNW KTSAGKEVIN KEDFVALERL  
TQGMDIQWMH VPGHSGFIGN EEADRLAREG AKQSED
```

## USAGE

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-0.2 mg/mL. Do not vortex.

## STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

## STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

## Image




Figure. SDS-PAGE

## **[IMPORTANT NOTE]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.