

Recombinant Caspase 7 (CASP7) Instruction Manual

SIPA782Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	23.2kDa
Accurate Molecular Mass	23kDa(Analysis of differences refer to the manual)
Residues & Tags	Ala24~Asp198 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	8.6
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
AKPDRSS IISSILLKKK RNASAGPVRT  
GRDRVPTYLY RMDFQKMGKC IIINNKNFDK ATGMDVRNGT DKDAGALFKC  
FQNLGFEVTV HNDSCAKMQ DLLRKASEED HSNSACFACV LLSHGEECLI  
YKGDGVTPIK DLTAHFRGDR CKTLLEKPKL FFIQACRGTE LDDGIQAD
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

G A T T C G C C A G C C A G A C C G C T C C T A T C A T C T C C T C T A T T C T C T T G A A G A A G A A G A A A T G C C T C T G C G G G C C C G T C A G G A C C G G C C G G A D C G A G T G C C C A C T T A T C T G T A C C G C A T G G A T T T C C A I




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.