

Recombinant Heat Shock Protein 40 (HSP40) Instruction Manual

SIPA900Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus, Cytoplasm
Predicted Molecular Mass	39.3kDa
Accurate Molecular Mass	39kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Ile340 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	8.8
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MGKDYQTLG LARGASDEEI KRAYRRQALR YHPDKNKEPG AEEKFKEIAE  
AYDVLSDPRK REIFDRYGEE GLKGSGPSGG SGGGANGTSF SYTFHGDPHA  
MFAEFFGGRN PFDTFQQRN GEEGMDIDDP FSGFPMGMGG FTNVNFGRSR  
SAQEPARKKQ DPPVTHDLRV SLEEIYSGCT KKMKISHKRL NPDGKSIRNE  
DKILTIEVKK GWKEGKITF PKEGDQTSNN IPADIVFVLK DKPHNIFKRD  
GSDVIYPARI SLREALCGCT VNVPTLDGRT IPVVFKDVIR PGMRRKVPGE  
GLPLPKTPEK RGDLIIEFEV IFPERIPQTS RTVLEQVLPI
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.


STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



MGEHTYDTQTGLARGAPDEEIKRATYRQALRYHPDANKKEFGAEEKFKETAEATDVLSDPRKREIFDRTGEEGLKGSQPSGGSGGGANGTSFSYTFHGDPHI





Figure. Gene Sequencing (Extract)

Image



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.