

Recombinant Protein Tyrosine Phosphatase Receptor Type C (CD45) Instruction Manual

SIPB015Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane
Predicted Molecular Mass	44.9kDa
Accurate Molecular Mass	44kDa(Analysis of differences refer to the manual)
Residues & Tags	Asp193~Lys575 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 90%
Isoelectric Point	6.4
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```

                                                                DAYLNASE
TTTTLSPSGSA VISTTTTIATT PSKPTCDEKY ANITVDYLYN KETKLFTAKL
NVNENVECGN NTCTNNEVHN LTECKNASVS ISHNSCTAPD KTLILDVPPG
VEKFQLHDCT QVEKADTTIC LKWKNIEFT CDTQNITYRF QCGNMIFDNK
EIKLENLEPE HEYKCDSEIL YNNHKFTNAS KIIKTDFGSP GEPQIIFCRS
EAAHQGVITW NPPQRSFHNF TLCYIKETEK DCLNLDKLI KYDLQNLKPY
TKYVLSLHAY IIAKVQRNGS AAMCHFTTKS APPSQVWMT VSMTSDNSMH
VKCRPPDRN GPHERYHLEV EAGNTLVRNE SHKNCDFRVK DLQYSTDYTF
KAYFHNGDYP GEPFILHHST SYNSK
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE


Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

GATGCTACCTTAAATGCTCTGAAACACCACTCTGAGGCTTCTGGAGGCTGCTTTTCAACACACAAATAGCTACTACTCCATCTAAGCCACATGTAAGAAAAATGTCAACATCCTGTGGATTCTTATATAA CAGGAACTAAATTATTTACAGCAAGCTA
D A Y L W A S E T T T L S P S G S A V I S T T T I A T T P S K P T C D E K Y A N I T V D Y L Y N K E T K L F T A K L



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.