

Recombinant Junctional Adhesion Molecule 1 (JAM1) Instruction Manual

SIPB174Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	n/a
Predicted Molecular Mass	23.9kDa
Accurate Molecular Mass	27kDa(Analysis of differences refer to the manual)
Residues & Tags	Phe64~Leu249 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	6.0
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
FDQGDTT RLVCYNMKIT ASYEDRVTFLL PTGITFKSVT  
REDTGTYTCM VSEEGGNSYG EVKVKLIVLV PPSKPTVNIP SSATIGNRAV  
LTCSEQDQSP PSEYTWFKDG IVMPTNPKST RAFSNSSYVL NPTTGELVFD  
PLSASDTGEY SCEARNGYGT PMTSNAVRME AVERN/GVIV AAVLVTLIL
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

