

Recombinant Granzyme K (GZMK) Instruction Manual

SIPB939Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Cytoplasm
Predicted Molecular Mass	24.1kDa
Accurate Molecular Mass	24kDa(Analysis of differences refer to the manual)
Residues & Tags	Gln44~Val227 with N-terminal His Tag
Buffer Formulation	100mMNaHCO ₃ , 500mMNaCl, pH8.3, containing 0.01% SKL, 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	9.6
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

CGGVLIHPQW VLTAAHCYSW FPRGHSPTVV LGAHSLSKNE PMKQTFEIKK QYRSKHI
FIPFSRLQSG SASHDIMLIK LRTAAELNKN VQLLHLGSKN YLRDGTKCQV
TGWGTTKPD L TASDTLREV TVTIISRKRC NSQSYYNHKP VITKDMICAG
DARGQKDSCK GDSGGPLICK GIFHALV

USAGE

Reconstitute in ddH₂O to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image


A horizontal bar chart representing a DNA or RNA sequence. The bar is composed of many small, multi-colored segments (red, green, blue, black) arranged in a regular, repeating pattern along the top of the page.

Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.