

Recombinant Myeloid Differentiation Factor 88 (MyD88) Instruction Manual

SIPB978Ga01

Chicken (Gallus)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	34.2kDa
Accurate Molecular Mass	34kDa(Analysis of differences refer to the manual)
Residues & Tags	Arg34~Pro299 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 90%
Isoelectric Point	5.9
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
                RLGLYLN PRAATAADWT  
ALAEKLGHDY LEIRRLEALP DPTAALLEEW QSRCPGGATV GQLELLRQL  
GRHDVLLLELG GSVEEDCKKY LRRKQQAEEQ PLQVPAVDSS VPKTSELMGI  
TTRDDPYGHG TEMFDAFICY CQKDLQFVQE MIRELEQTEF KKLKCVFDRD  
VLPGTCVWSI SGELIERRCR RMVVVISDDY LESDECDFQT KFALSLSPGA  
RLKRLIPVKC KTMKNEFPSI LRFITICDYT NPCTKKWFWT RLAKSLLLP
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

