

# Recombinant Activating Transcription Factor 3 (ATF3) Instruction Manual

## SIPC111Mu01

### Mus musculus (Mouse)

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Nucleus
<b>Predicted Molecular Mass</b>	50.7kDa
<b>Accurate Molecular Mass</b>	51kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Met1~Ser181 with N-terminal His and GST Tag
<b>Buffer Formulation</b>	100mMNaHCO <sub>3</sub> , 500mMNaCl, pH8.3, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 95%
<b>Isoelectric Point</b>	8.8
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

## SEQUENCE

```
MMLQHPGQVS ASEVSATAIV PCLSPPGSLV FEDFANLTPF VKEELRFAIQ  
NKHLCHRMSS ALESVTVNNR PLEMSVTKSE AAPEEDERKR RRRERNKIAA  
AKCRNKKKEK TECLQKESEK LESVNAELKA QIEELKNEKQ HLIYMLNLHR  
PTCIVRAQNG RTPEDERNLF IQQIKEGTLQ S
```

## USAGE

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## STORAGE

