

Recombinant Thyroid Hormone Receptor Alpha (THRa) Instruction Manual

SIPC412Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus
Predicted Molecular Mass	58.5kDa
Accurate Molecular Mass	59kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Pro490 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01%skl, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	6.4
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MEQKPSKVEC GSDPEENSAR SPDGKRKRKN GQCSLKTSMS GYIPSYLDKD
EQCVVCGDKA TGYHYRCITC EGCKGFFRRT IQKNLHPTYS CKYDSCCVID
KITRNQCQLC RFKKCIAVGM AMDLVLDDSK RVAKRKLIEQ NRERRRKEEM
IRSLQQRPEP TPEEWDLIHI ATEAHRSTNA QGSHWKQRRK FLPDDIGQSP
IVSMPDGDKV DLEAFSEFTK IITPAITRVV DFAKKLPMFS ELPCEQIIL
LKGCCMEIMS LRAAVRYDPE SDTLTSLGEM AVKREQLKNG GLGVVSDAIF
ELGKSLSAFN LDDTEVALLQ AVLLMSTDRS GLLCVDKIEK SQEAYLLAFE
HYVNRKHNI PHFWPKLLMK EREVQSSILY KGAAAEGRPG GSLGVHPEGQ
QLLGMHVVGQ PQVRQLEQQL GEAGSLQGPV LQHSPKSPQ QRLLELLHRS
GILHARAVCG EDDSSEADSP SSSEEEPEVC EDLAGNAASP
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.