

Recombinant Cytochrome P450 1B1 (CYP1B1) Instruction Manual

SIPD110Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	17.9kDa
Accurate Molecular Mass	18/16kDa(Analysis of differences refer to the manual)
Residues & Tags	Asp374~Phe516 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	7.9
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
      DQPNLPY VLAFLYEAMR FSSFVPVTIP  
HATTANTSVL GYHIPKDTV  FVNQWSVNHD PLKWPNPENF DPARFLDKDG  
LINKDLTSRV MIFSVGKRRC IGEELSKMQL FLFISILAHQ CDFRANPNEP  
AKMNFSGYGLT IKPKSF
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

G A A T T C G A C C A G C C C A A C C T G C C C T A T G T C C T G G C C T T C C T T A T G A A G C C A T G C G C T T C C A G C T T T G T G C C T G T C A C T A T T C C T C A T G C C A C C A C T G C C A A C A C C T C T G T C T T G G G C T A C C A C A T T C




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.