

Recombinant Diazepam Binding Inhibitor (DBI) Instruction Manual

SIPD264Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Endoplasmic reticulum lumen
Predicted Molecular Mass	11.1kDa
Accurate Molecular Mass	13kDa(Analysis of differences refer to the manual)
Residues & Tags	Ser2~Ile87 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	8.0
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

SQAEFDKAA EEVKRLKTQP TDEEMLFIYS HFKQATVGDV NTDRPGLLDL
KGGKAKWDSWN KLKGTSKESA MKTYVEKVDE LKKKYGI

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE


Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

TCTGAGCTGATTGACAAAGGCGTGAAGGAGGAGGAGGCTAAGACTCAGCCAGTGTGAGAGAGTCTGTTCACTCAGTCACTCAACAGCTAAGTGGGCGATGATACAGATGGGCGGGCTTTGGACTCAGGGCAAGCCAGTGGGACTGTGGACAACTGAAAGG
S Q A E F D K A A E E V K R L K T Q P T D E E W L F I T S R F X Q A T V G D V N T D R P G L I D L K G K A K V D S W N X L E G



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.