

Recombinant Interleukin 17C (IL17C) Instruction Manual

SIPD343Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	21.0kDa
Accurate Molecular Mass	25kDa(Analysis of differences refer to the manual)
Residues & Tags	His19~Val197 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	8.2
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
HH DPSLRGHPHS HGTPHCYSAE ELPLGQAPPH  
LLARGAKWGQ ALPVALVSSL EAASHRGRHE RPSATTQCPV LRPEEVLEAD  
THQRSISPWR YRVDTDEDRY PQLAFAECL CRGCIDARTG RETAALNSVR  
LLQSLLVLRG RPCSRDGSGL PTPGAFAFHT EFIHVPVGCT CVLPRSV
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.