

Recombinant Troponin T Type 1, Slow Skeletal (TNNT1) Instruction Manual

SIPD410Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	34.2kDa
Accurate Molecular Mass	40kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Gly259 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 90%
Isoelectric Point	5.7
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MSDTEEQEYE EEQPEEEAAE EEEEEPEEPE PVAEPEEERP KPSRPVVPPL  
IPPKIPEGER VDFDDIHRKR MEKDLELQTL LIDVHFEQRK KEEELVALK  
ERIERRRSER AEQQRFRTEK ERERQAKLAE EKMRKEEEEA KKRAEDDAKK  
KKVLSNMGAAH FGGYLVKAEQ KRGKRQTGRE MKVRILSERK KPLDIDYMGE  
EQLREKAQEL SDWIHQLESE KFDLMAKLLQ QKYEINVLVN RISHAQKFRK  
GAGKGRVGG
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

