

Recombinant Myosin VA (MYO5A) Instruction Manual

SIPD772Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Lysosome, Cytoplasm, Golgi apparatus, Endoplasmic reticulum lumen, Endosome
Predicted Molecular Mass	40.7kDa
Accurate Molecular Mass	41kDa(Analysis of differences refer to the manual)
Residues & Tags	Ser1529~Val1853 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	8.3
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
SL LTSTINSIKK VLKKGDDFE  
TVSFWLSNTC RFLHCLKQYS GEEGFMKHNT SRQNEHCLTN FDLAEYRQVL  
SDLAIQIYQQ LVRVLENILQ PMIVSGMLEH ETIQGVSGVK PTGLRKRTSS  
IADEGTYTLD SILRQLNSFH SVMCQHGM DP ELIKQVVKQM FYIVGAILTN  
NLLLRRKDMCS WSKGMQIRYN VSQLEEWLRD KNLMNSGAKE TLEPLIQAAQ  
LLQVKKKTDD DAEAICSMCN ALTTAQIVKV LNLYTPVNEF EERVSVSFIR  
TIQMRLRDRK DSPQLLMDAK HIFPVTFPFN PSSLALETIQ IPASLGLGFI  
ARV
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.