

Recombinant Histone Deacetylase 6 (HDAC6) Instruction Manual

SIPE313Hu03

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus, Cytoplasm
Predicted Molecular Mass	42.7kDa
Accurate Molecular Mass	43kDa(Analysis of differences refer to the manual)
Residues & Tags	Ser479~Arg835 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 80%
Isoelectric Point	6.6
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
SR TGLVYDQNMN NHCNLWDSHH  
PEVPQRILRI MCRLEELGLA GRCLTLTPRP ATEAELLTCH SAEYVGHRLA  
TEKMKTRELH RESSNFDSIY ICPSTFACAQ LATGAACRLV EAVLSGEVLN  
GAAVVRPPGH HAEQDAACGF CFFNSVAVAA RHAQTISGHA LRILIVDWDV  
HHGNGTQHMF EDDPSVLYVS LHRYDHGTFE PMGDEGASSQ IGRAAGTGFT  
VNVAWNGPRM GDADYLAAWH RLVLPYAYEF NPELVLVSAQ FDAARGDPLG  
GCQVSPEGYA HLTHLLMGLA SGRIILILEG GYNLTSISES MAACTRSLLG  
DPPPLLLTPR PPLSGALASI TETIQVHRRY WRSLR
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.