

Recombinant Sirtuin 1 (SIRT1)

Instruction Manual

SIPE319Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus
Predicted Molecular Mass	30.0kDa
Accurate Molecular Mass	30kDa(Analysis of differences refer to the manual)
Residues & Tags	Glu236~Ala490 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	5.9
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
                                EDAVK LLQECKKIIV
LTGAGVSVSC GIPDFRSRDG IYARLAVDFP DLPDPQAMFD IEYFRKDPRP
FFKFAKEIYP GQFQPSLCHK FIALSDKEGK LLRNYTQNIID TLEQVAGIQR
ILQCHGSFAT ASCLICKYKV DCEAVRGDIF NQVVPRCPRC PADEPLAIMK
PEIVFFGENL PEQFHRAMKY DKDEVDLLIV IGSSLKVRPV ALIPSSIPHE
VPQILINREP LPHLHFDVEL LGDCDVIINE LCHRLGGEYA
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.

STORAGE

