

# Recombinant CREB Binding Protein (CREBBP) Instruction Manual

SIPE684Ra01

*Rattus norvegicus* (Rat)

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Nucleus, Cytoplasm
<b>Predicted Molecular Mass</b>	25.4kDa
<b>Accurate Molecular Mass</b>	26kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Gln2216~Leu2442 with N-terminal His Tag
<b>Buffer Formulation</b>	20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 95%
<b>Isoelectric Point</b>	7.0
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

## SEQUENCE

```
          QQSSA  SLAGGMAGHS  QFQQPQGP GG  YAPAMQQQRM
QQHLPIQGSS  MGQMAAPMGQ  LGQMGQPGLG  ADSTPNIQQA  LQQRILQQQQ
MKQQIGSPGQ  PNPMSPPQQHM  LSGQPQASHL  PGQQIATSLS  NQVRSPAPVQ
SPRPQSQPPH  SSPSPRIQPQ  PSPHHVSPQT  GSPHPGLAVT  MASSMDQGHL
GNPEQSAMPLP  QLNTPNRSAL  SSELSLVGDT  TGD TLEKFVE  GL
```

## USAGE

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-1.0 mg/mL. Do not vortex.


## STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

## STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

### Image



SDS-PAGE Image

### **[IMPORTANT NOTE]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.