

Recombinant Tight Junction Protein 2 (TJP2) Instruction Manual

SIPF496Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane, Nucleus
Predicted Molecular Mass	44.2KDa
Accurate Molecular Mass	49kDa(Analysis of differences refer to the manual)
Residues & Tags	Gly307~Ala669 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	5.1
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
GVLL MKSRANEEYG LRLGSQIFVK EMTRTGLATK DGNLHEGDII
L KINGTVTEN MSLTDARKLI EKSRGKLQLV VLRDSQQLI NIPSLNDS
EIEDISEIES NRSFSPEERR HQYSDYDYHS SSEKLKERPS SREDTPSRLS
RMGATPTPFK STGDIAGTVV PETNKEPRYQ EDPPAPQPKA APRTFLRPSP
EDEAIYGPNT KMVRFKKGDS VGLRLAGGND VGIFVAGIQE GTSAEQEG LQ
EGDQILKVNT QDFRGLVRED AVLYLLEIPK GEMVTILAQS RADVYRDILA
CGRGDSFFIR SHFECEKETP QSLAFTRGEV FRVVD TLYDG KLG NWLAVRI
GNELEKGLIP NKSRAEQMA
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.


STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



15% SDS-PAGE

SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.