

Recombinant Utrophin (UTRN) Instruction Manual

SIPF595Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Membrane, Cytoplasm
Predicted Molecular Mass	32.5kDa
Accurate Molecular Mass	33kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Glu252 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 97%
Isoelectric Point	6.8
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MAKYGEHEAS PDNGQNEFSD IIKRSDEHN DVQKKTFTKW INARFSKSGK  
PPINDMFTDL KDGRKLLDLL EGLTGTS LPK ERGSTRVHAL NNVNRVLQVL  
HQNVELVNI GGTDIVDGNH KLTGLLLWSI ILHWQVKDVM KDVMSDLQQT  
NSEKILLSWV RQTRTPYSQV NVLNFTTSWT DGLAFNAVLH RHKPD LFSWD  
KVVKMSP IER LEHAFSKAQT YLGIEKLLDP EDVAVQLPDK KSIIMYLTSL  
FE
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.