

# Recombinant Microfibrillar Associated Protein 4 (MFAP4) Instruction Manual

**SIPF819Hu01**

**Homo sapiens (Human)**

|                                 |                                                                            |
|---------------------------------|----------------------------------------------------------------------------|
| <b>Source</b>                   | Prokaryotic expression                                                     |
| <b>Host</b>                     | E.coli                                                                     |
| <b>Endotoxin Level</b>          | <1.0EU per 1µg (determined by the LAL method)                              |
| <b>Subcellular Location</b>     | n/a                                                                        |
| <b>Predicted Molecular Mass</b> | 30.2kDa                                                                    |
| <b>Accurate Molecular Mass</b>  | 32kDa(Analysis of differences refer to the manual)                         |
| <b>Residues &amp; Tags</b>      | Val22~Ala255 with                                                          |
| <b>Buffer Formulation</b>       | PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5%<br>Trehalose and Proclin300. |
| <b>Traits</b>                   | Freeze-dried powder                                                        |
| <b>Purity</b>                   | > 95%                                                                      |
| <b>Isoelectric Point</b>        | n/a                                                                        |
| <b>Applications</b>             | Positive Control; Immunogen; SDS-PAGE; WB.                                 |

## **USAGE**

Reconstitute in PBS or others.

## **STORAGE**


Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

## **STABILITY**

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no

obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**Image**



**[IMPORTANT NOTE]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.