

Recombinant Keratin 25 (KRT25) Instruction Manual

SIPS706Mu01

Mus musculus (Mouse)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Cytoplasm
Predicted Molecular Mass	52.6KDa
Accurate Molecular Mass	53kDa(Analysis of differences refer to the manual)
Residues & Tags	Met1~Asn446 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	4.7
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
MSLRLSSGSR RSYARPSTGS LRGASFGAGN ACGVAGIGSG FSCAFGGSSST  
GGNTGVANSC AGFTVNEGGL LSGNEKVTMQ NLNDRLASYL DNVQALQEAN  
ADLEQKIKGW YEKFGPGSCR GLDHDYSRYF PIIDDLKNQI ITSTTSNANA  
VLQIDNARLT ADDFRLKYEN ELALHQSVEA DVNGLRRVLD EITLCRTDLE  
IQYETLSEEL TYLKKNHKEE MQALQCAAGG NVNVEMNAAP GVDLTVLLNN  
MRAEYEALAE QNRRDAEAWF QEKASLQQQ ITEDVGATTS ARNELTEMKR  
TLQTLIEIQ SLLATKHSLE CSLTETEGNY CTQLAQIQAQ ISALEEQLHQ  
VRTETEGQKL EYEQLLNKVA HLEKEIETYC LLIGGDEGAC KSSSYKSKDY  
GSGNAGNQIK DPVKAIVVKK VLEEVDQRSK ILTTRLHSLE EKSQSN
```

USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.
Do not vortex.


STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.

